



velvet XP white



**Comfort and hygiene for
everyday clinical care**



Hospital



Dental



Primary
care



Long-term
care



Veterinary

Sempermed is a HARPS Global Brand.*

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet XP white

sempermed® velvet XP white is a powder-free nitrile examination glove designed for reliable barrier protection and precise handling in clinical environments. The glove combines a comfortable fit with dependable tactile sensitivity and a clean white appearance for professional use. An additional variant with extended cuff length (velvet XP white+) provides increased coverage for situations requiring enhanced splash protection.

- White color**
 For a professional and hygienic appearance.
- Textured fingertips**
 Offer secure grip at all times, even in damp conditions.
- Elastic nitrile formulation**
 For comfortable fit and tactile control.
- Extended cuff option (velvet XP white+)**
 For increased coverage and additional splash protection.

| Technical data | | | | | | | |
|---|--|--------|------|--------------|--------------|-----------------|-----------------|
| Type | non-sterile, powder-free examination glove for single use | | | | | | |
| Color | white | | | | | | |
| Size / Overall length as per EN 455-2 | XS, S, M, L, XL <i>standard: median ≥ 240 mm</i> <i>velvet XP+ white: median ≥ 265 mm</i> | | | | | | |
| Wall thickness measured in single layer | <table border="0"> <tr> <td>Finger</td> <td>Palm</td> </tr> <tr> <td>min. 0.07 mm</td> <td>min. 0.06 mm</td> </tr> <tr> <td>typical 0.10 mm</td> <td>typical 0.07 mm</td> </tr> </table> | Finger | Palm | min. 0.07 mm | min. 0.06 mm | typical 0.10 mm | typical 0.07 mm |
| Finger | Palm | | | | | | |
| min. 0.07 mm | min. 0.06 mm | | | | | | |
| typical 0.10 mm | typical 0.07 mm | | | | | | |
| Barrier performance | AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement) | | | | | | |
| Material | Nitrile Butadiene Rubber (NBR) not made with natural rubber latex | | | | | | |
| Glove surface | finger textured | | | | | | |
| Glove shape | ambidextrous with rolled rim | | | | | | |
| Force at break as per EN 455-2 | median ≥ 6 N | | | | | | |
| Shelf life | 5 years in original packaging if stored according to product specification | | | | | | |
| Labeling & conformity to standards | EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M | | | | | | |

| Purpose | |
|--|--|
| Single-use medical device class I as set out in MD Regulation (EU) 2017/745 | |
| Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425 | |
| Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials | |

| Packaging | |
|-------------------|--|
| Dimensions | 150 Dispenser: 240L x 122W x 65H mm Carton: 340L x 252W x 250H mm |
| | 200 Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm |
| Available | in dispenser box XS, S, M, L 150 / XL 135 pcs in dispenser box XS, S, M, L 200 / XL 180 pcs 10 dispenser boxes per shipping carton |

| Article numbers | | | | |
|-----------------|-----------------|------------|------------------|------------|
| | velvet XP white | | velvet XP+ white | |
| XS | 200 pcs | 3000016021 | 150 pcs | 3000016532 |
| S | 200 pcs | 3000016022 | 150 pcs | 3000016533 |
| M | 200 pcs | 3000016023 | 150 pcs | 3000016534 |
| L | 200 pcs | 3000016024 | 150 pcs | 3000016535 |
| XL | 180 pcs | 3000016025 | 135 pcs | 3000016536 |

*Sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd., an affiliate of HARPS Global Pte. Ltd.
 **Please find below chemical resistance tables for our different glove brands: <https://www.sempermed.com/chemical-resistance-brands>
IMPORTANT NOTES: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **HARPS Global is not liable for incorrect use of the gloves nor will HARPS Global be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classifications set out in this datasheet are up to date as of the date of publication. Where typical technical details are listed, they are lot-specific values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. ©Copyright 2026 HARPS Global