

S **semper**
guard

sapphire™ XP black



Sleek and Professional with Enhanced Performance



Medical



Beauty &
Tattoo



Safety &
Protection



Automotive



Janitorial
& Sanitation



Food &
Beverage

*semperguard is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





sapphire™ XP black

Semperguard® sapphire™ XP black exam nitrile gloves are designed to offer great grip and tactile sensitivity, making them ideal for a variety of tasks that require precision and control. The sleek black color delivers a professional appearance while effectively concealing color transfer, stains, inks and other messes, maintaining a cleaner and more polished appearance, throughout use. Ideal for industries such as beauty, laboratory, food services, janitorial-sanitation, and safety, where aesthetics and cleanliness are a priority.

- Sleek and professional appearance.**
 Hides stains caused from grease, dirt, dyes, and ink
- Added versatility**
 Exam grade with extended sizing
- Tested against 36 chemotherapy and hazardous drugs, fentanyl, and gastric acid**
 Meets USP 800 requirements for PPE compliance for chemotherapy
- Safe for food contact applications**
 Meets FDA 21 CFR Part 177 requirements for food handling

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Black	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.9 mils	Palm 2.8 mils
	based on average production value measured in single layer	
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	575% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, ASTM D6978	

Packaging

DIMENSIONS

Dispenser: 125L x 245W x 42H mm
Carton: 223L x 260W x 255H mm

Reorder numbers / Glove Counts*

	GSBNF101	X-Small	100 gloves/bx	10 bxs/cs
	GSBNF102	Small	100 gloves/bx	10 bxs/cs
	GSBNF103	Medium	100 gloves/bx	10 bxs/cs
	GSBNF104	Large	100 gloves/bx	10 bxs/cs
	GSBNF105	X-Large	100 gloves/bx	10 bxs/cs
	GSBNF106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks. For USP <800> chemotherapy use, change gloves every 20-30 minutes.

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global