

S semper
guard

vinyl classic™



**Simply Safe,
and Simply Classic**



**General
Purpose**



**Janitorial &
Sanitation**



**Safety &
Protection**



**Beauty &
Tattoo**



**Food &
Beverage**

*semperguard is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





vinyl classic™

Semperguard® vinyl classic™ powdered gloves are designed to provide reliable performance for everyday tasks. These gloves are suitable for users in the safety and janitorial-sanitation industry who require a practical and dependable solution. Semperguard® vinyl classic™ offers functionality and value, featuring a smooth surface and a formulation free from DEHP and DINP.

- Compliant with California Prop 65**

Formulated without DEHP and DINP plasticizers.

- Powdered for easy donning**

Helping to reduce friction and improving consistent fit in demanding environments.

Technical data





TYPE	Non-sterile vinyl examination and disposable protective glove	
GRADE	General Purpose	
COLOR	Clear	
MATERIAL	Polyvinyl Chloride (PVC)	
GLOVE SURFACE	Smooth	
COUNTRY OF ORIGIN	China	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.9mils	Palm 2.6mils
	based on average production value measured in single layer	
TENSILE STRENGTH <i>(before aging)</i>	11MPa (min.)	
ELONGATION <i>(before aging)</i>	300% (min.)	
BARRIER PERFORMANCE	AQL < 4.0 (pinholes)	
TESTING & CONFORMITY TO STANDARDS	EN ISO 21420	

Packaging

DIMENSIONS

Dispenser: 200L x 110W x 60H mm
Carton: 310L x 225W x 210H mm

Reorder numbers / Glove Counts*

	VP102	Small	100 gloves/bx	10 bxs/cs
	VP103	Medium	100 gloves/bx	10 bxs/cs
	VP104	Large	100 gloves/bx	10 bxs/cs
	VP105	X-Large	100 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing vinyl gloves between tasks.

PEOPLE, PROTECTING PEOPLE

rev 06/2026 US

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global