

S **semper**
med

latex edition IC™



Comfort and Strength



Ambulatory
Surgical
Centers



Dental



Veterinary



Long Term
Care



Primary
Care



Non-medical
Lab

*sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





latex edition IC™

Sempermed® latex edition IC™ exam gloves are designed to provide the perfect balance of protection, tactile sensitivity and value. The advanced manufacturing technology and synthetic polymer inner coating creates a glove that is strong, flexible and goes on with ease.

- Specially formulated with **comfort and protection** in mind.
- Synthetic polymer inner coating allows for **easy donning**.
- The textured surface finish offers a **strong grip** in wet or dry environments.

Technical data






TYPE	Non-sterile latex examination and disposable protective glove	
GRADE	Medical	
COLOR	Natural	
MATERIAL	Made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 94mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 5.5mils	Palm 4.7mils
based on average production value measured in single layer		
TENSILE STRENGTH <i>(before aging)</i>	20MPa (avg.)	
ELONGATION <i>(before aging)</i>	709% (avg.)	
BARRIER PERFORMANCE	AQL < 1.5 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D3578, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5	

Packaging

DIMENSIONS

Dispenser: 240L x 122W x 65H mm
Carton: 340L x 252W x 250H mm

Reorder numbers / Glove Counts*

	SCLT101	X-Small	100 gloves/bx	10 bxs/cs
	SCLT102	Small	100 gloves/bx	10 bxs/cs
	SCLT103	Medium	100 gloves/bx	10 bxs/cs
	SCLT104	Large	100 gloves/bx	10 bxs/cs
	SCLT105	X-Large	100 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

 Licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gogloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global