

semper  
med

rose™



## Care, Comfort and a Touch of Color



Hospitals



Dental



Veterinary



Ambulatory  
Surgical  
Centers



Primary  
Care



Beauty &  
Tattoo

sempermed is a registered trademark of  
HARPS Investment Asia Pte. Ltd.

CONTACT US! [gloveinfo@harpsglobal.com](mailto:gloveinfo@harpsglobal.com) | [www.gotgloves.com](http://www.gotgloves.com)

HARPS USA, Inc.  
3001 Executive Drive, Suite 330  
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.  
1 Harbourfront Avenue, #14-07 Keppel Bay Tower  
Singapore 098632, Singapore





# rose™

Sempermed® rose™ exam nitrile gloves are designed to elevate comfort and protection in any setting, from exam rooms to tattoo shops to day spas. These gloves feature a soothing aloe coating that helps prevent skin irritation, keeping hands moisturized and comfortable during extended use. The unique formulation of these gloves enhances protection while maintaining exceptional tactile sensitivity, allowing for precise handling without compromising comfort.

- Aloe & vitamin E inner coating**  
 Helps moisturize and care for dry hands, supporting skin comfort during extended glove use.
- Fun unique color**  
 Distinctive pink color creates a friendly, professional appearance and helps bring a positive touch to clinical and wellness environments.
- Protects users and patients from the transmission of infections**  
 Designed to provide protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

## Technical data






<b>TYPE</b>	Non-sterile nitrile examination and disposable protective glove	
<b>GRADE</b>	Medical	
<b>COLOR</b>	Pink	
<b>MATERIAL</b>	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
<b>GLOVE SURFACE</b>	Textured	
<b>COUNTRY OF ORIGIN</b>	Malaysia	
<b>SHelf LIFE</b>	5 years in original packaging if stored according to product specification	
<b>SIZE</b>	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
<b>WALL THICKNESS</b>	Finger 3.9mils	Palm 2.8mils
	based on average production value measured in single layer	
<b>TENSILE STRENGTH</b> <i>(before aging)</i>	31MPa (avg.)	
<b>ELONGATION</b> <i>(before aging)</i>	542% (avg.)	
<b>BARRIER PERFORMANCE</b>	AQL < 1.5 (pinholes)	
<b>SOCIAL COMPLIANCE</b>	amfori BSCI, EcoVadis	
<b>TESTING &amp; CONFORMITY TO STANDARDS</b>	ASTM D6319, EN ISO 21420	

## Packaging

### DIMENSIONS

Dispenser: 245L x 125W x 74H mm  
Carton: 385L x 255W x 255H mm

## Reorder numbers / Glove Counts\*

	<b>SMNR201</b>	X-Small	200 gloves/bx	10 bxs/cs
	<b>SMNR202</b>	Small	200 gloves/bx	10 bxs/cs
	<b>SMNR203</b>	Medium	200 gloves/bx	10 bxs/cs
	<b>SMNR204</b>	Large	200 gloves/bx	10 bxs/cs
	<b>SMNR205</b>	X-Large	180 gloves/bx	10 bxs/cs

## Certified Claims



SAMPLES: Available by request - ask your sales rep!

\*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

## PEOPLE, PROTECTING PEOPLE

**IMPORTANT NOTE:** The latest product information is at [www.gotgloves.com/products/](http://www.gotgloves.com/products/). Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global