

semper
med

velvet™



Soft as Velvet
and Light as a Feather



Hospitals



Dental



Veterinary



Long Term
Care



Primary
Care



Ambulatory
Surgical
Centers

sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





velvet™

Sempermed® velvet™ is a lightweight, soft examination and protective glove. It is a facility-wide glove that is suitable for a wide range of applications. This glove is the perfect latex-free alternative for sensitive skin. The soft nitrile formulation allows the glove to conform to the hand like a second skin, combining optimal fit and comfort with unrestricted hand motion.

● **Soft as velvet**

The soft formulation of the glove offers enhanced comfort and fit.

● **Your daily choice**

Suitable for a wide variety of medical applications where protection is critical.

● **Tested against 56 chemotherapy and hazardous drugs, fentanyl, and gastric acid**

● **Protects users and patients from the transmission of infections**

Designed to provide a high level barrier of protection against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Violet Blue	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 3.3mils	Palm 2.4mils based on average production value measured in single layer
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	575% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, ASTM D6978	

Packaging

DIMENSIONS	300ct	Dispenser: 245L x 128W x 83H mm Carton: 430L x 266W x 257H mm
	100ct	Dispenser: 245L x 125W x 42H mm Carton: 260L x 223W x 255H mm

Reorder numbers / Glove Counts*

SMNP301 SMNS101	X-Small	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP302 SMNS102	Small	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP303 SMNS103	Medium	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP304 SMNS104	Large	300 gloves/bx 100 gloves/bx	10 bxs/cs
SMNP305 SMNS105	X-Large	250 gloves/bx 100 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

Velvet 100ct. is licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global