



# vinyl comfort™



## Comfort Without Compromise



Long Term Care



Veterinary



Lab & Pharma



Primary Care

sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd.

CONTACT US! [gloveinfo@harpsglobal.com](mailto:gloveinfo@harpsglobal.com) | [www.gotgloves.com](http://www.gotgloves.com)

HARPS USA, Inc.  
3001 Executive Drive, Suite 330  
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.  
1 Harbourfront Avenue, #14-07 Keppel Bay Tower  
Singapore 098632, Singapore





# vinyl comfort™

Sempermed® vinyl comfort™ exam gloves are a reliable and cost-effective solution for professionals seeking a latex-free alternative for routine applications. These single-use vinyl gloves are manufactured without DEHP and DINP plasticizers and are suitable for a wide range of applications.

- Compliant with California Proposition 65**  
 Formulated without DEHP and DINP plasticizers.
- Polymer coated**  
 Helping to reduce friction and improving consistent fit in demanding environments.

## Technical data

<b>TYPE</b>	Non-sterile vinyl examination and disposable protective glove	
<b>GRADE</b>	Medical	
<b>COLOR</b>	Clear	
<b>MATERIAL</b>	Polyvinyl Chloride (PVC)	
<b>GLOVE SURFACE</b>	Smooth	
<b>COUNTRY OF ORIGIN</b>	China	
<b>SHELF LIFE</b>	5 years in original packaging if stored according to product specification	
<b>SIZE</b>	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
<b>WALL THICKNESS</b>	Finger 3.3mils	Palm 3.1mils based on average production value measured in single layer
<b>TENSILE STRENGTH</b> <i>(before aging)</i>	18MPa (avg.)	
<b>ELONGATION</b> <i>(before aging)</i>	406% (avg.)	
<b>BARRIER PERFORMANCE</b>	AQL < 1.5 (pinholes)	
<b>TESTING &amp; CONFORMITY TO STANDARDS</b>	ASTM D5250, EN ISO 21420	

## Packaging

<b>DIMENSIONS</b>	Dispenser: 240L x 122W x 60H mm Carton: 315L x 252W x 256H mm
-------------------	--

## Reorder numbers / Glove Counts\*

	<b>EVNP102</b>	Small	100 gloves/bx	10 bxs/cs
	<b>EVNP103</b>	Medium	100 gloves/bx	10 bxs/cs
	<b>EVNP104</b>	Large	100 gloves/bx	10 bxs/cs
	<b>EVNP105</b>	X-Large	100 gloves/bx	10 bxs/cs

## Certified Claims



SAMPLES: Available by request - ask your sales rep!  
\*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing vinyl gloves between tasks.

Licensed for use in Canada

## PEOPLE, PROTECTING PEOPLE

**IMPORTANT NOTE:** The latest product information is at [www.gotgloves.com/products/](http://www.gotgloves.com/products/). Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global