



vinyl syn-stretch



The Smarter Stretch



Long Term
Care



Veterinary



Lab &
Pharma



Primary
Care

sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





vinyl syn-stretch

Sempermed® vinyl syn-stretch gloves are designed for users seeking enhanced comfort and flexibility compared to conventional vinyl gloves. With improved elasticity and tear resistance, these single-use gloves provide a soft feel and comfortable fit for a wide range of routine applications.

Compliant with California Proposition 65

Formulated without DEHP and DINP plasticizers.

Polymer coated

Helping to reduce friction and improving consistent fit in demanding environments.

Technical data

TYPE	Non-sterile vinyl examination and disposable protective glove	
GRADE	Medical	
COLOR	Yellow	
MATERIAL	Polyvinyl Chloride (PVC)	
GLOVE SURFACE	Smooth	
COUNTRY OF ORIGIN	China	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 230mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 4.1mils Palm 3.3mils	based on average production value measured in single layer
TENSILE STRENGTH (before aging)	23MPa (avg.)	
ELONGATION (before aging)	469% (avg.)	
BARRIER PERFORMANCE	AQL < 1.5 (pinholes)	
TESTING & CONFORMITY TO STANDARDS	ASTM D5250, EN ISO 21420	

Packaging

DIMENSIONS

Dispenser: 240L x 122W x 60H mm
Carton: 315L x 252W x 256H mm

Reorder numbers / Glove Counts*

SCVNP102	Small	100 gloves/bx	10 bxs/cs
SCVNP103	Medium	100 gloves/bx	10 bxs/cs
SCVNP104	Large	100 gloves/bx	10 bxs/cs
SCVNP105	X-Large	100 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing vinyl gloves between tasks.



Licensed for use in Canada

PEOPLE, PROTECTING PEOPLE

rev 06/2026 US

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global