

semper
med

safe™



**Strong and
Durable**



Hospitals



Lab &
Pharma



Veterinary



First
Responder



Primary
Care



Beauty &
Tattoo

*sempermed is a registered trademark of
HARPS Investment Asia Pte. Ltd.*

CONTACT US! gloveinfo@harpsglobal.com | www.gotgloves.com

HARPS USA, Inc.
3001 Executive Drive, Suite 330
Clearwater, FL 33762

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Singapore





safe™

Sempermed® safe™ exam nitrile gloves combine durability, comfort, and reliable protection for professional use. The nitrile material offers resistance to a broad range of industrial chemicals, making these gloves suitable for demanding applications across multiple industries.

- Good chemical resistance**
 Against a wide array of chemicals, offers the user enhanced protection and safety.
- Tested against 56 chemotherapy and hazardous drugs, fentanyl, and gastric acid**
- The textured surface on the glove's fingertips offers a strong grip in wet or dry environments.
- Protects users and patients from the transmission of infections**
 Designed to provide a high level of protection and barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination and infection in a clinical setting.
- ESD tested**
 According to ANSI/ESD STM 15.1

Technical data

TYPE	Non-sterile nitrile examination and disposable protective glove	
GRADE	Medical	
COLOR	Blue	
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
GLOVE SURFACE	Textured	
COUNTRY OF ORIGIN	Malaysia	
SHELF LIFE	5 years in original packaging if stored according to product specification	
SIZE	Length: 240mm (min) Width: 95mm +/- 5mm (based on size M)	
WALL THICKNESS	Finger 5.1mils	Palm 3.5mils
	based on average production value measured in single layer	
TENSILE STRENGTH <i>(before aging)</i>	31MPa (avg.)	
ELONGATION <i>(before aging)</i>	564% (avg.)	
BARRIER PERFORMANCE	AQL < 1.0 (pinholes)	
SOCIAL COMPLIANCE	amfori BSCI, RGA, WRAP, EcoVadis	
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, EN ISO 374-1, EN ISO 374-5, ISO 18889, ASTM D6978	

Packaging

DIMENSIONS

Dispenser: 245L x 123W x 54H mm
Carton: 283L x 257W x 261H mm

Reorder numbers / Glove Counts*

	BLNF102	Small	100 gloves/bx	10 bxs/cs
	BLNF103	Medium	100 gloves/bx	10 bxs/cs
	BLNF104	Large	100 gloves/bx	10 bxs/cs
	BLNF105	X-Large	100 gloves/bx	10 bxs/cs
	BLNF106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



avg. 10° D under worst-case test conditions.

SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"

Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks.

PEOPLE, PROTECTING PEOPLE

rev 06/2026 US

IMPORTANT NOTE: The latest product information is at www.gotgloves.com/products/. Failure to observe the information provided, in particular with regard to chemical resistance, frequency of use, and tolerability of the gloves, can result in personal injury and/or material damage. HARPS Global is not liable for any incorrect use of the gloves, improper storage and handling, or allergic reactions. In case of damage or unusual signs of wear and tear, gloves should not be used. If there are any doubts, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors, and amendments. ©2026 HARPS Global