



skySense black



**Everyday protection in
professional black**



**Safety &
Protection**



**Cleaning &
Sanitation**



**Food &
Beverage**



**Beauty &
Body-Art**



Medical

*Semperguard is a HARPS Global brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





skySense black

Semperguard® skysense black examination and protective gloves are a lightweight and reliable solution for hygiene-sensitive tasks requiring precision and control. The thin construction supports good tactile sensitivity and wearing comfort, while the finger-textured surface enhances grip in wet and dry conditions.

The black color provides a professional appearance and helps conceal stains, making the glove particularly suitable for food handling, catering, cleaning, and beauty or body-art applications. The classification as medical device supports use in hygiene-critical environments where barrier protection is required.

Black color for professional use

Provides a clean, professional appearance while helping to conceal stains

Finger-textured for improved grip

Enhances control in wet and dry conditions for precise handling.

Barrier protection for hygienic applications

Glove material is resistant to blood and other body fluids and has been tested against bacteria, fungi and viruses, supporting hygienic procedures.

Comfortable and flexible

Thin construction supports tactile sensitivity and wearing comfort during extended use.

Technical data

Type	non-sterile, powder-free, examination and disposable protective glove	
Color	black	
Size / Overall length as per EN 455-2	XS, S, M, L, XL, / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.05 mm typical 0.07 mm	Palm min. 0.05 mm typical 0.06 mm
Barrier performance as per EN 455-1	AQL ≤ 1.5	
Material	nitrile butadiene rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	flat with rolled rim, ambidextrous	
Force at break as per EN 455-2	median ≥ 6 N	
Shelf life	3 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III (protection against chemical substances for limited time)** as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 245L x 125W x 42H mm
	Carton: 260L x 223W x 255H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs 10 dispenser boxes per shipping carton

Article numbers

XS	100 pcs	3000016469
S	100 pcs	3000016470
M	100 pcs	3000016471
L	100 pcs	3000016472
XL	90 pcs	3000016473