


semperclean MC



**Anatomical comfort
for demanding tasks**



**Safety &
Protection**



**Lab &
Pharma**



**Primary
care**



Hospital

*Semperguard is a HARPS Global brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +6567976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





semperclean MC

Semperclean® MC is a high-quality natural latex examination glove developed for demanding tasks in laboratory, pharmaceutical, and healthcare environments. Its anatomical shape, based on surgical glove formers, supports a natural fit and high wearing comfort, while the extended cuff provides additional coverage compared to standard examination gloves.

The QuickDon™ inner coating enables smooth and easy donning, supporting efficient glove changes in fast-paced settings such as laboratory work or hospital preparation procedures. The micro-roughened surface enhances grip and handling precision, making the glove well suited for detailed tasks and controlled working environments, including non-sterile medical applications such as cleaning, disinfection, and preparation prior to sterile procedures.

■ Anatomical shape for high wearing comfort

Based on surgical glove formers for a natural fit and reduced hand fatigue.

■ Extended cuff for additional coverage

Longer than a standard examination glove for enhanced wrist coverage.

■ Easy donning with QuickDon™ inner coating

Supports fast and smooth glove changes in high-frequency tasks.

■ Micro-roughened surface for precise grip

Supports secure handling and control during detailed tasks.

Technical data

Type	non-sterile, powder-free, examination and disposable protective glove		
Color	natural white		
Size / Overall length as per EN 455-2	6.0 / 6.5 7.0 / 7.5 / 8.0 8.5 / 9.0	/	median ≥ 260 mm median ≥ 270 mm median ≥ 280 mm
Wall thickness measured in single layer	Palm Finger	0.17 – 0.23 mm max. 0.27 mm	
Barrier performance as per EN 455-1	AQL ≤ 1.5		
Material	natural rubber (NR)		
Glove surface	micro-rough		
Glove shape	fully anatomical with curved fingers and rolled rim		
Force at break as per EN 455-2	typical ≥ 12 Newton		
Shelf life	3 years in original packaging if stored according to product specification		
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 421, ISO 2859-1, ASTM F1671/F1671M		

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 280 L x 110 W x 185 H mm Carton: 283 L x 382 W x 570 H mm
Available	50 pairs in dispenser box 500 pairs in transport carton

Article numbers

size 6	814720602
size 6½	814720622
size 7	814720702
size 7½	814720722
size 8	814720802
size 8½	814720822
size 9	814720902