

semper
guard

iconiX



Visible Protection Against Invisible Risks



First
Responder



Safety &
Protection



ASC



Lab &
Pharma



Pest



Hospital

*Semperguard is a HARPS Global brand.**

CONTACT US!
gloveinfo@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +6567976560

HARPS USA, Inc.
3001 Executive Drive, Suite 330,
Clearwater, FL 33762
Tel. +1 800 366 9545





iconiX

Semperguard® iconiX examination gloves are engineered for professionals who demand uncompromising protection. These durable nitrile gloves deliver dependable barrier performance for clinical, laboratory, pharmaceutical, and emergency response applications. The two-tone design supports the visual detection of breaches and helps prevent accidental reuse, while the extended cuff provides enhanced forearm coverage. Tested against 56 chemotherapy and hazardous drugs, including Carmustine and Thiotepe, Semperguard® iconiX gloves provide trusted protection when handling highly hazardous substances. They also meet ISO 18889 requirements for use with pesticides and are NFPA 1950(1999)-2025 certified to enhance safety for emergency medical operations.

- **Tested for 56 Chemotherapy and hazardous drugs,** according to ASTM D6978, including Thiotepe (>240 min) and Carmustine (77.6 min)
- **NPFA 1950(1999)-2025 certified** for use as an emergency medical examination glove
- **Meets G1 protection according to ISO 18889** suitable for handling diluted pesticides
- **Extended cuff** for added protection and for additional wrist and forearm coverage
- **Special damp-don inner coating** facilitates faster donning, even with damp hands

Technical data

TYPE	Non-sterile, powder-free, examination and disposable protective glove
GRADE	Exam
COLOR	Dual color blue
MATERIAL	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex
GLOVE SURFACE	Finger textured
COUNTRY OF ORIGIN	Malaysia
SHELF LIFE	3 years in original packaging if stored according to product specification
SIZE	Length: 295mm (avg.) Width: median 95 mm ± 5 mm (based on size M)
WALL THICKNESS	Palm 5.1 mils Finger 7.1 mils based on average production value measured in single layer
TENSILE STRENGTH (after aging)	36 MPa (avg.)
ELONGATION (after aging)	528% (avg.)
BARRIER PERFORMANCE	AQL ≤ 1.0 (pinholes)
QUALITY CERTIFICATION	Responsible Glove Alliance (RGA), amfori BSCI, WRAP, EcoVadis
TESTING & CONFORMITY TO STANDARDS	ASTM D6319, EN 455, EN ISO 21420, ISO 11193-1, EN ISO 374-1, EN ISO 374-5, ISO 18889, ASTM D6978, NFPA 1950 (1999)-2025

Packaging

DIMENSIONS	Dispenser: 270 L x 130 W x 65 H mm Carton: 336 L x 268 W x 280 H mm
------------	--

Reorder numbers / Glove Counts*

ICNX102	Small	100 gloves/bx	10 bxs/cs
ICNX103	Medium	100 gloves/bx	10 bxs/cs
ICNX104	Large	100 gloves/bx	10 bxs/cs
ICNX105	X-Large	90 gloves/bx	10 bxs/cs
ICNX106	XX-Large	90 gloves/bx	10 bxs/cs

Certified Claims



SAMPLES: Available by request - ask your sales rep!

*glove counts "by weight"
Gloves are intended for single use. For maximum protection we recommend changing nitrile gloves between tasks. For USP <800> chemotherapy use, change gloves every 20–30 minutes.