

S semper
guard

iconix



Visible Protection Against Invisible Risks



**Safety &
Protection**



**Lab &
Pharma**



**Primary
care**



Hospital

*Semperguard is a HARPS Global brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +6567976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





iconiX

Semperguard® iconiX gloves combine protection and functionality and have been tested against 56 chemotherapy and hazardous drugs. The light-blue exterior allows short notes to be written directly on the glove, supporting quick identification and documentation during use. The dual-color design helps distinguish the outside from the inside of the glove and reduces the risk of unintentional re-use. The gloves provide chemical resistance across an extended list of chemicals and feature a longer cuff for additional coverage. Certified as both a Medical Device (MD) and Personal Protective Equipment (PPE), semperguard® iconiX gloves are designed for use in demanding environments.

- **Tested against 56 chemotherapy and hazardous drugs**
for protection against a broad range of challenging substances.**
- **Designed with a longer cuff**
to provide added coverage and extended wrist protection.
- **Special damp-don inner coating**
facilitates faster donning, even with damp hands.
- **The innovative dual-color design**
provides a clear visual cue to prevent unintentional re-use.
- **Light-blue outside**
allows short notes to be written directly on the glove for quick identification and documentation.

Technical data

Type	non-sterile, powder-free examination and disposable protective glove		
Color	dual-color blue		
Size / Overall length as per EN 455-2	S, M, L, XL, XXL / median 295 mm ± 5 mm		
Wall thickness measured in single layer	Finger min. typical	0.15 mm 0.16 mm	Palm 0.12 mm 0.13 mm
Barrier performance as per EN 455-1	AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement)		
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex		
Glove surface	finger textured		
Glove shape	flat with rolled rim, ambidextrous		
Force at break as per EN 455-2	typical ≥ 13 N		
Shelf life	3 years in original packaging if stored according to product specification		
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 421, ISO 2859-1, ASTM F1671/F1671M		

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 245 L x 125 W x 42 H mm Carton: 260 L x 223 W x 255 H mm
Available	in dispenser box S, M, L 50 / XL, XXL 45 pcs 10 dispenser boxes per shipping carton

Article numbers

S	50 pcs	3000016541
M	50 pcs	3000016542
L	50 pcs	3000016543
XL	45 pcs	3000016544
XXL	45 pcs	3000016545