

S semper
med

velvet XP



**Fits like a
second skin**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global brand**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet XP

sempermed® velvet XP is a powder-free nitrile examination glove designed for everyday use across a wide range of applications. It combines a soft and flexible material with reliable barrier protection and tactile control, supporting comfortable wear that fits like a second skin during routine procedures. Compared to the standard velvet version, it offers slightly increased material strength for users who prefer an added level of protection in daily practice.

■ Soft and flexible material

For a comfortable fit and tactile control.

■ Xtra Protection (XP)

Slightly increased material strength for users preferring an added level of protection.

■ Textured fingertips

For reliable grip, even in damp conditions.

■ Protects users and patients from the transmission of infections**

Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	violet blue	
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.07 mm typical 0.10 mm	Palm min. 0.06 mm typical 0.07 mm
Barrier performance	AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	flat with rolled rim	
Force at break as per EN 455-2	typical ≥ 7 N	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)** as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions	100 pcs	Dispenser: 245L x 125W x 42H mm Carton: 260L x 223W x 255H mm
	200 pcs	Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs in dispenser box XS, S, M, L 200 / XL 180 pcs 10 dispenser boxes per shipping carton	

Article numbers

XS	100 pcs	3000015971	200 pcs	3000015991
S	100 pcs	3000015972	200 pcs	3000015992
M	100 pcs	3000015973	200 pcs	3000015993
L	100 pcs	3000015974	200 pcs	3000015994
XL	90 pcs	3000015975	180 pcs	3000015995