



velvet sterile



Sterile & soft as velvet

Sempermed is a HARPS Global brand.*



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet sterile

With sempermed® velvet sterile, the sempermed brand offers a sterile version of its lightweight and exceptionally soft examination and protective glove. Designed for applications where enhanced hygiene standards are required, it is a reliable choice for a wide range of non-surgical procedures.

As the ideal latex-free alternative, sempermed® velvet sterile combines outstanding wearing comfort with reliable mechanical strength. Its soft nitrile formulation provides a second-skin feel and is particularly suitable for users with sensitive skin.

The sterile design supports stricter infection prevention protocols and offers a dependable alternative to surgical gloves for non-surgical procedures. This ensures a high level of safety and confidence in environments where minimizing the risk of contamination is essential.

■ Soft as velvet

The soft material results in high wearing comfort and a comfortable fit.

■ Your daily choice

Suitable for various application fields as per the intended use.

■ Protects users and patients from the transmission of infections.

Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	sterile, powder-free examination glove for single use	
Color	violet blue	
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.05 mm typical 0.09 mm	Palm min. 0.05 mm typical 0.06 mm
Barrier performance	AQL ≤ 0.65 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	ambidextrous with rolled rim	
Force at break as per EN 455-2	median ≥ 6 N	
Shelf life	3 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, EN 556-1 ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class Is as set out in MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III (protection against chemical substances for limited time)** as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 155L x 135W x 250H mm Carton: 480L x 285W x 261H mm
Available	in dispenser box XS, S, M, L, XL 50 pairs 6 dispenser boxes per shipping carton

Article numbers

XS	50 pairs	3000016234
S	50 pairs	3000016235
M	50 pairs	3000016236
L	50 pairs	3000016237
XL	50 pairs	3000016238