

velvet



**Soft as velvet
and light as a feather**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632, Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet

sempermed® velvet is a powder-free nitrile examination glove designed for everyday use across a wide range of applications. It combines a soft and flexible material with reliable barrier protection and tactile control, supporting comfortable wear during routine procedures. The glove is suitable for users who require consistent performance in daily clinical and care environments.

■ Soft as velvet

The soft material results in high wearing comfort and a comfortable fit.

■ Your daily choice

Suitable for various application fields as per the intended use.

■ Protects users and patients from the transmission of infections.

Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	violet blue	
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.05 mm typical 0.09 mm	Palm min. 0.05 mm typical 0.06 mm
Barrier performance	AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	ambidextrous with rolled rim	
Force at break as per EN 455-2	median ≥ 6 N	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	100 pcs	Dispenser: 245L x 125W x 42H mm Carton: 260L x 223W x 255H mm
	200 pcs	Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs in dispenser box XS, S, M, L 200 / XL 180 pcs 10 dispenser boxes per shipping carton	

Article numbers

XS	100 pcs	3000015965	200 pcs	3000015946
S	100 pcs	3000015966	200 pcs	3000015947
M	100 pcs	3000015967	200 pcs	3000015948
L	100 pcs	3000015968	200 pcs	3000015949
XL	90 pcs	3000015969	180 pcs	3000015950