



latex comfort



LATEX ALL-ROUNDER



Hospital



Dental



Primary
care



Long-term
care



Veterinary

Sempermed is a HARPS Global Brand.*

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





latex comfort

sempermed® latex comfort is a powder-free latex examination glove intended for single use in medical and occupational settings. With its lightweight design, it offers good tactile control and wearing comfort while meeting the required standards for barrier protection. The glove is suitable for routine procedures and high-frequency use.

High tactile sensitivity

Due to the elastic properties of natural latex.

Lightweight design

For comfortable wear and efficient use.

Textured fingertips

For secure grip in wet and dry conditions.

Protects users and patients from the transmission of infections

Glove material is resistant to blood and other body fluids, and protects against bacteria, fungi and viruses including SARS-CoV-2, Hepatitis B & C and HIV.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	natural white	
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median ≥ 240 mm	
Wall thickness measured in single layer	Finger min. 0.09 mm typical 0.11 mm	Palm min. 0.08 mm typical 0.10 mm
Barrier performance	AQL ≤ 1.5	
Material	Natural Rubber (NR)	
Glove surface	finger textured	
Glove shape	ambidextrous with rolled rim	
Force at break as per EN 455-2	median ≥ 6 N	
Shelf life	3 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type C), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D3578 (except stress at 500% elongation), ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 240L x 122W x 65H mm Carton: 340L x 252W x 250H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs 10 dispenser boxes per shipping carton

Article numbers

XS	100 pcs	3000016036
S	100 pcs	3000016037
M	100 pcs	3000016038
L	100 pcs	3000016039
XL	90 pcs	3000016040