

S semper
med

safe



Strong and durable



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





safe

Feel safe with sempermed® safe, a powder-free nitrile examination glove designed for reliable barrier protection across a wide range of applications. Its durable construction supports everyday use in medical and industrial environments, combining mechanical strength with a comfortable fit and dependable grip. The glove is suitable for users requiring consistent performance and broad application versatility.

■ The textured surface at the glove's fingertips

Offers strong grip in wet or dry conditions.

■ The good chemical resistance**

Against a wide array of industrial chemicals offers the user enhanced protection and safety.

■ Strong and durable nitrile material

Makes the glove suitable for many different uses.

■ The force at break

Is significantly higher than the EN 455-2 standard requirement, offering the user an enhanced protection.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	ocean blue	
Size / Overall length as per EN 455-2	S, M, L, XL	/ median ≥ 240 mm
Wall thickness measured in single layer	Finger min. 0.10 mm typical 0.13 mm	Palm min. 0.08 mm typical 0.09 mm
Barrier performance	AQL 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	flat with rolled rim	
Force at break as per EN 455-2	typical ≥ 9 N	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, EN 421, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 245L x 123W x 54H mm Carton: 280L x 254W x 255H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs 10 dispenser boxes per shipping carton

Article numbers

S	100 pcs	3000015996
M	100 pcs	3000015997
L	100 pcs	3000015998
XL	90 pcs	3000015999