

S semper
med

safe+



Extra long. Extra safe.

*Sempermed is a HARPS Global Brand.**



Hospital



Primary
care



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





safe+

The sempermed® safe+ provides reliable hand protection when extra barrier and splash resistance are needed. Thanks to the glove's length, wrists and forearms are better protected, which makes it a suitable choice for glove users who work with risk patients or handle substances that can be harmful to the skin.

Enhanced safety

Additional splash protection thanks to the extra-long cuff.

Chemical resistance

Barrier protection against certain liquids and chemicals.

Work with confidence

Finger texturing and good tactile sensitivity.

Suitable for sensitive skin

Nitrile material helps to prevent type-1 allergies.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	blue	
Size / Overall length as per EN 455-2	XS,S, M, L, XL / median $\geq 295 \text{ mm} \pm 5 \text{ mm}$	
Wall thickness measured in single layer	Finger min. 0.12 mm typical 0.14 mm	Palm min. 0.08 mm typical 0.09 mm
Barrier performance	AQL ≤ 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	flat with rolled rim	
Force at break as per EN 455-2	typical $\geq 10 \text{ N}$	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 421, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	Dispenser: 280L x 125W x 60H mm Carton: 310L x 260W x 90H mm
Available	in dispenser box XS, S, M, L 100 / XL 90 pcs 10 dispenser boxes per shipping carton

Article numbers

XS	100 pcs	3000016026
S	100 pcs	3000016027
M	100 pcs	3000016028
L	100 pcs	3000016029
XL	90 pcs	3000016030