

S semper
med

skySense



Your daily choice of nitrile examination gloves

In line with **Regulation (EU) 2022/1031** on the International Procurement Instrument – IPI, skySense can be sourced from Malaysia and China under the same article code.



Hospital



Dental



Primary
care



Long-term
care



Veterinary

Sempermed is a HARPS Global brand.*

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348

sempermed® skySense offers a lightweight, comfortable and resilient solution for various applications.

The lightweight nitrile examination and protective gloves are available in three colors: blue, violet blue and white. The gloves are not only strong, but also designed to support unrestricted movement and ease of use. They deliver functionality and comfort at great value.

- **Your daily choice:** suitable for a wide variety of medical applications where protection is essential.
- **Tested for protection against chemotherapy drugs.****
- **Protects users and patients from the transmission of infections.** Designed to provide a reliable barrier against contaminants, bacteria, and viruses, reducing the risk of cross-contamination in clinical settings.

Technical data

Type	non-sterile, powder-free examination glove for single use
Color	blue / violet blue / white
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median 240 mm
Barrier performance as per EN 455-1	AQL 1.5
Material	Nitrile Butadiene Rubber (NBR)
Glove surface	finger textured
Glove shape	flat with rolled rim, ambidextrous
Force at break as per EN 455-2	median ≥ 6 N
Shelf life in original packaging if stored according to product specification	3 years
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III (protection against chemical substances for limited time)** as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

available in 100 and 200 pieces per dispenser box
in dispenser box XS, S, M, L 100 / XL 90 pieces
in dispenser box XS, S, M, L 200 / XL 180 pieces
10 dispenser boxes per shipping carton

Dimensions & Loading Qtys.

DIMENSIONS	Dispenser 245L x 125W x 42H mm 100 pcs. Carton 260L x 223W x 255H mm
200 pcs.	Dispenser 245L x 125W x 74H mm Carton 380L x 258W x 255H mm
LOADING QUANTITIES	20ft 2,033 / 1,187 cartons 40ft 4,175 / 2,483 cartons 100/200 pcs. 40ft HC 4,639 / 2,760 cartons

DM1/17 Version 3 08/2026 EN

Article Numbers	skySense	skySense VB	skySense white
XS 100pc 3000016142	200pc 3000015792	100pc 3000016463 200pc 3000015800	200pc 3000016489
S 100pc 3000016143	200pc 3000015793	100pc 3000016464 200pc 3000015821	200pc 3000016490
M 100pc 3000016144	200pc 3000015794	100pc 3000016465 200pc 3000015822	200pc 3000016491
L 100pc 3000016145	200pc 3000015795	100pc 3000016466 200pc 3000015823	200pc 3000016492
XL 90pc 3000016146	180pc 3000015796	90pc 3000016467 180pc 3000015824	180pc 3000016493

*Sempermed is a registered trademark of HARPS Investment Asia Pte. Ltd., an affiliate of HARPS Global Pte. Ltd.

**Please find the chemical resistance tables for our different glove brands at: <https://www.harpsglobal.com/chemical-resistance-brands>

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **HARPS Global is not liable for incorrect use of the gloves nor will HARPS Global be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **HARPS Global bears no liability for allergic reactions.** *Copyright 2026 HARPS Global