



velvet XP white



**Comfort and hygiene for
everyday clinical care**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

*Sempermed is a HARPS Global Brand.**

CONTACT US!
sempermed@harpsglobal.com
www.harpsglobal.com

HARPS Investment Asia Pte. Ltd.
1 Harbourfront Avenue, #14-07 Keppel Bay Tower
Singapore 098632
Tel.: +65 67976560

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348





velvet XP white

sempermed® velvet XP white is a powder-free nitrile examination glove designed for reliable barrier protection and precise handling in clinical environments. The glove combines a comfortable fit with dependable tactile sensitivity and a clean white appearance for professional use. An additional variant with extended cuff length (velvet XP+ white) provides increased coverage for situations requiring enhanced splash protection.

White color

For a professional and hygienic appearance.

Textured fingertips

Offer secure grip at all times, even in damp conditions.

Elastic nitrile formulation

For comfortable fit and tactile control.

Extended cuff option (velvet XP+ white)

For increased coverage and additional splash protection.

Technical data

Type	non-sterile, powder-free examination glove for single use	
Color	white	
Size / Overall length as per EN 455-2	XS, S, M, L, XL	standard: median $\geq$ 240 mm velvet XP+ white: median $\geq$ 265 mm
Wall thickness measured in single layer	Finger min. 0.07 mm typical 0.10 mm	Palm min. 0.06 mm typical 0.07 mm
Barrier performance	AQL $\leq$ 1.0 (exceeding the EN 455-1 barrier requirement)	
Material	Nitrile Butadiene Rubber (NBR) not made with natural rubber latex	
Glove surface	finger textured	
Glove shape	ambidextrous with rolled rim	
Force at break as per EN 455-2	median $\geq$ 6 N	
Shelf life	5 years in original packaging if stored according to product specification	
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M	

Purpose

Single-use medical device class I as set out in
MD Regulation (EU) 2017/745



Single-use protective glove – PPE category III
(protection against chemical substances for limited time)**
as set out in PPE Regulation (EU) 2016/425



Suitable for food contact as set out in Regulation (EC)
1935/2004 on Food Contact Materials



Packaging

Dimensions	150	Dispenser: 240L x 122W x 65H mm Carton: 340L x 252W x 250H mm
	200	Dispenser: 245L x 125W x 74H mm Carton: 380L x 258W x 255H mm
Available	in dispenser box XS, S, M, L 150 / XL 135 pcs in dispenser box XS, S, M, L 200 / XL 180 pcs 10 dispenser boxes per shipping carton	

Article numbers

	velvet XP white		velvet XP+ white	
XS	200 pcs	3000016021	150 pcs	3000016532
S	200 pcs	3000016022	150 pcs	3000016533
M	200 pcs	3000016023	150 pcs	3000016534
L	200 pcs	3000016024	150 pcs	3000016535
XL	180 pcs	3000016025	135 pcs	3000016536